




# How to Use a Metered-Dose Inhaler with an Aerochamber (Spacer)




1. Shake the medicine.



2. Insert the mouthpiece of the inhaler into the rubber-sealed end of the Aerochamber.



3. Breathe all the air out of your lungs then put the Aerochamber into your mouth between your teeth. Make a tight seal around the mouthpiece with your lips.



4. Press the metered-dose inhaler down once to release a spray of medicine. The medicine will be trapped in the spacer. Breathe in slowly and deeply.



5. Hold your breath for at least 5 to 10 seconds. Breathe out slowly.

Copyright ©2001 McKesson Health Solutions LLC. All rights reserved.